



30K Burs*

FOR THE STRAIGHTSHOT® M5 MICRODEBRIDER

*Maximum Speed
for ENT Surgery*



*Pending clearance

30K Burs*

FOR THE STRAIGHTSHOT® M5 MICRODEBRIDER

Maximum Speed for ENT Surgery.

Finding the right bur for the job has never been easier. With the largest selection of burs for frontal sinus surgery, you can choose a bur for the specific anatomic needs of each patient or procedure.



Ordering Info

Product #	Description	Length
1884015RTD	30K Tapered Diamond Choanal Atresia Bur, 15°	13 cm
1883040BLD	30K Bullet Diamond Bur, 3.0 mm, 40°	13 cm
1883040BLC	30K Bullet Cutting Bur, 3.0 mm, 40°	13 cm
1884040RTC	30K Reverse Tapered Cutting Bur, 4.0 mm, 40°	13 cm
1883040BRC	30K Frontal Finesse Barrel Cutting Bur, 3.0 mm, 40°	13 cm
1883655BRC	30K Barrel Cutting Bur, 3.6 mm, 55°	13 cm
1884070RTD	30K Reverse Tapered Diamond Bur, 4.0 mm, 70°	13 cm
1884070RTC	30K Reverse Tapered Cutting Bur, 4.0 mm, 70°	13 cm
1883070BLD	30K Bullet Diamond Bur, 3.0 mm, 70°	13 cm
1883070BLC	30K Bullet Cutting Bur, 3.0 mm, 70°	13 cm

30K Burs* have an entirely new design and multiple technology innovations, offering:

- **Speeds up to 30,000 rpm**
- **Best selection available for frontal sinus drilling**
- **Integrated irrigation**
- **Distal suction**

Rx only. Refer to product instruction manual/package insert for instructions, warnings, precautions, and contraindications.

*Pending clearance

For further information, please call Medtronic ENT at 800.874.5797 or consult Medtronic's website at www.medtronic.com or www.medtronic.com/30Kburs.

Medtronic, Inc. ENT

6743 Southpoint Drive N
Jacksonville, FL 32216
USA
Toll free: (800) 874-5797
Telephone: (904) 296-9600
Fax: (800) 678-3995

International Telephone Numbers

Adriatic Region 385-1-488-1120
Argentina 54-11-4898-5700
Australia 1-800-668-670
Baltic Region 37-1-67560226
Belgium 32-2456-09-09
Brazil 55-11-2182-9200
Canada 1-800-217-1617
Chile 56-2-2655-5110
China 86-21-50800998
Colombia 57-1-742-7300
Czech Republic 420-2-9657-9580
France 33-470-679-800
Germany 49-2159-8149-209
Greece 30-210-67-79-099
Hong Kong 852-2919-1312
Hungary 36-188-90600
India 91-22-26836733
Israel 972-9-972-4400
Italy 39-02-24137-324
Japan 81-3-6430-2017
Korea 82-2-3404-3600
Lebanon 961-1-370-670
Luxembourg 32-2456-09-09
Malaysia 60-37-953-4800
Mexico 52-55-11-02-90-30
Netherlands 31-45-566-8800
Poland 48-22-4656900
Russian Fed. 7-495-580-73-77
Singapore 65-6776-6255
South Africa 27-11-466-1820
Spain 34-91-625-05-40
Taiwan 886-2-2183-6000
Thailand 662-232-7400
UK 44-1923-205-166



Medtronic