

Integra[®] MicroFrance[®]

Vaughan Suction Sinus Instruments

Developed in conjunction with Winston Vaughan, MD



These instruments have been designed to evacuate blood from the surgical site to provide visibility in the sinuses.

- Suction mushroom punches are designed with the suction at the tip of the instrument
- A cleaning tool allows soft tissue to be removed during surgery
- Suction picks help aid in the dissection of closed cells while evacuating blood



Ordering Information

Reference	Description
MCEN110	Punch Cleaner, 4 and 6mm, Double-Ended, Working Length 121mm (4¾")
MCEN111	Suction Pick, 3mm, 110° Forward, Working Length 106mm (4⅛")
MCEN112	Suction Pick, 3mm, 90° Back, Working Length 105mm (4⅛")
MCEN114	Suction Cup Forceps, 2.8mm 30° Up Angled, Working Length 123mm (4⅞")
MCEN115-6	Malleable Suction, 6Fr, 125mm (4⅞")
MCEN115-8	Malleable Suction, 8Fr, 120mm (4¾")
MCEN115-9	Malleable Suction, 9Fr, 125mm (4⅞")
MCEN116	Suction Mushroom Punch, 4mm, Straight, Working Length 117mm (4⅝")
MCEN117	Suction Mushroom Punch, 4mm, 60° Up, Working Length 125mm (4⅞")
MCEN118	Suction Mushroom Punch, 6mm, Straight, Working Length 117mm (4⅝")
MCEN120	Suction Sickle Knife, 3mm, 105mm (4⅛")
3717002	Instrument Tray, 52cm L x 25cm W x 2.8cm D, (20½" L x 9¾" W x 1⅛" D)



The MicroFrance® brand is a premier line of ENT instruments, designed with world-renowned ENT surgeons to facilitate new techniques and improve patient outcomes.

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

■ Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.

■ Warning: Applicable laws restrict these products to sale by or on the order of a physician.

■ Consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions, and instructions for use.

For more information or to place an order, please contact:

United States, Canada, Asia, Pacific, Latin America

USA 800-654-2873 ■ 888-980-7742 fax

International +1 609-936-5400 ■ +1 609-750-4259 fax

integralife.com/contact



Integra MicroFrance SAS
Le Pavillon, 03160
Saint-Aubin-Le-Monial ■ FRANCE