

Integra[®] MicroFrance[®]

Myringotomy Instruments


A Comprehensive Set



THE COMPLETE SOLUTION FOR

Myringotomy Procedures

This comprehensive set of instruments provides everything you need for myringotomy procedures. It includes our most popular specula, blades, forceps, and other instruments to accommodate your patient needs.



1 Myringotomy Instrument Set

Reference: 3754010*

Set includes the following:

Reference	Description	Qty
3721026	House Mini Forceps, 4.5mm Delicate Serrated Jaws	1
3724013	Buck Ear Curette, Angled, Dull, Size 1, 2.0mm	1
3724055	Billeau Ear Loop, Small, 2.0mm Loop	1
3722058	Barbara Pick, Straight Point, Malleable Shaft	1

Ebony Specula

Reference	Description	Qty
3724401	Ear Speculum, Oval, Oblique, Size 1, 3 X 4mm, Ebony	1
3724403	Ear Speculum, Oval, Oblique, Size 3, 4 X 5mm, Ebony	1
3724405	Ear Speculum, Oval, Oblique, Size 5, 5 X 6mm, Ebony	1

Suction Tubes


Reference	Description	Qty
3724172	Baron Suction Tube, 75mm Working Length, 3Fr	1
3724173	Baron Suction Tube, 75mm Working Length, 5Fr	1
3724175	Baron Suction Tube, 75mm Working Length, 7Fr	1

Tray

Reference	Description	Qty
MCEN136SC	Instrument Tray, 285mm L x 175mm W x 50mm D	1



3724401 3724403 3724405
3 x 4mm 4 x 5mm 5 x 6mm




MCEN136SC

*This set is only available in the US. If you are not in the US, please order the above items individually.


INTEGRA[®]
LIMIT UNCERTAINTY



3721026



3724013



3724055



3722058



3 FR: 3724172
5 FR: 3724173
7 FR: 3724175



2

Select a Handle and Blade

Myringotomy Blade Handles†

Reference	Description	Qty
1352321	E/Z Chuck Blade Handle	1
1352325	LiteTouch Uni-Handle-Black Hardco at Aluminum	1

An alternative to the E/Z Chuck Blade Handle


Single-Use Myringotomy Blades†

Reference	Description	Qty
1458230	Micro Otologic Blade, Upcutting, Angled	5
1458231	Micro Otologic Blade, Downcutting, Angled	5
1458232	Micro Otologic Blade, Razor, Straight	5
1458233	Micro Otologic Blade, Spear, Angled	5
1458234	Micro Otologic Blade, Spear, Straight	5
1458275	Myringotomy Blade, Narrow, Angled	5
1458276	Myringotomy Blade, Spear Tip, Angled	5
1411000	Myringotomy Knife, Spear Blade, Integral Handle	10
1411050	Myringotomy Knife, Lance Blade, Integral Handle	10

†Use color-coded blue and yellow dots to match compatible handles and blades.



1352321



1352325




1458230

1458233



1458231



1458234

1458232



1458276

1458275




1411000

1411050

3 Add Disposables

Microlite Disposable Suction Needles

Reference	Description	Qty
1359254	Microlite Disposable Suction Needle, 14G, Red	10
1359256	Microlite Disposable Suction Needle, 16G, Yellow	10
1359258	Microlite Disposable Suction Needle, 18G, Black	10
1359260	Microlite Disposable Suction Needle, 20G, Green	10
1359262	Microlite Disposable Suction Needle, 22G, Clear	10
1359264	Microlite Disposable Suction Needle, 22G, Blue	10
1358270	Disposable Suction Handle with Bypass Control	10



Disposable Myringotomy Kit

Reference	Description	Qty
0208	This sterile, single use kit is a convenient all-in-one solution including: • 5 sizes of ear specula • 2 sizes of curettes • 1 myringotomy knife • 1 suction needle 16 gauge	10



4


Include any Optional Accessories

Instruments

Reference	Description	Qty
3714010	Lathbury Applicator, Fine Twist 20mm Tip, 1.2mm Diam	1
3722132	Straight Tip Needle, Angled Shaft	1
3724028	Day Ear Hook, Small, 3.0mm, Hook Angled 90°	1
3755000	Backhaus Towel Forceps, 88mm (3.5")	1
3755001	Backhaus Towel Forceps, 133mm (5.2")	1

Curettes

Reference	Description	Qty
3724017	Buck Ear Curette, Straight, Dull, Size 0, 1.5mm	1
3724018	Buck Ear Curette, Straight, Dull, Size 1, 2.0mm	1



Forceps

Reference	Description	Qty
3711002	Lucae Dressing Forceps, Bayonet Style	1
3721025	Alligator Forceps, Very Fine, 6.0mm Serrated Jaws	1
3721028	Alligator Forceps, Mini, 4.5mm Smooth Jaws	1
3721036	Alligator Forceps, Very Fine, 4.0mm Serrated Jaws	1



Specula & Holders

Reference	Description	Qty
3724093	Boucheron Speculum, Round, Size 1, 3.0mm	1
3724094	Boucheron Speculum, Round, Size 2, 3.5mm	1
3724095	Boucheron Speculum, Round, Size 3, 4.5mm	1
3724096	Boucheron Speculum, Round, Size 4, 5.0mm	1
3724316	Shea Speculum Holder	1
1362279	Universal Speculum Holder Set**	1



3724093
3.0mm

3724094
3.5mm

3724095
4.5mm

3724096
5.0mm



** 1362279 has the following components that must be ordered separately outside the US: 1362275 Universal Flexible Arm, 1352282 Table Bracket Blade, and 1352277 Quik-Lock Speculum Mount.

The MicroFrance® brand is a premier line of ENT instruments, designed with world-renowned ENT surgeons to facilitate new techniques and improve patient outcomes.

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Always refer to the appropriate instructions for use for complete clinical instructions.
- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.

For more information or to place an order, please contact:

United States, Canada, Asia, Pacific, Latin America

USA 800-654-2873 ▪ 888-980-7742 fax

International +1 609-936-5400 ▪ +1 609-750-4259 fax

integralife.com/contact



Integra MicroFrance SAS
Le Pavillon, 03160
Saint-Aubin-le-Monial ▪ France



Integra York PA, Inc.
589 Davies Drive
York, PA 17402 ▪ USA