

Integra[®] Jarit[®]

Video Assisted Thoracoscopic Surgery



Limit Uncertainty

with the Brands you Trust.



Table of Contents

Clamps4

Forceps8

Needle Holders16

Scissors18

Dissector19

Node Graspers.....20

Suction21

Knot Tier/Pushers24

Index.....25

Guarantee/Warranty

All Integra® Jarit® Instruments are unconditionally guaranteed against defects in material and workmanship. Integra Jarit will, at its option and without charge, either repair or replace any instrument which Integra Jarit determines to be defective in material or workmanship when used for its intended surgical purpose.

All instruments sold by Integra Jarit are warranted only to the original purchaser and only against defects in workmanship or materials which under the intended use render such instruments inoperable.

All Integra Jarit Carb-Bite™ and Integra Jarit CARB-EDGE™ instruments are further warranted for a period of three (3) years under normal use. Should any adjustment become necessary during this warranty period, CARB-EDGE Scissors will be resharpened, and Carb-Bite Needle Holder jaw inserts will be replaced without charge.

All Integra Jarit Supercut™ Scissors are further warranted for a period of one year under normal use. Should any resharpening become necessary during this warranty period, Integra Jarit will provide this service without charge. A nominal charge will be made for sharpening after expiration of this warranty. Integra Jarit Supercut Scissors will be date-stamped for control.

This Guarantee/Warranty shall be void in the event any Integra Jarit instrument is used for purposes other than as originally intended, is not maintained or cleaned properly, is damaged as a result of misuse or accident or is repaired or altered by unauthorized persons.

The above is a limited warranty, and it is the only warranty made by Integra Jarit. Integra Jarit makes no other warranty, express or implied, including any warranty of merchantability or fitness for a particular purpose, as well as any express or implied warranty to patients. Integra Jarit shall have no liability for consequential, exemplary or incidental damages even if it is aware of the possibility of such damages. No warranty or guarantee may be created by any act or statement nor may this warranty be modified in any way, except as a result of a writing signed by an officer of Jarit. These limitations on the creation or modification of this warranty may not be waived or modified orally or by any conduct. The stated express warranty is in lieu of all liabilities or obligations of Integra Jarit arising out of, or in connection with, the delivery, use or performance of any Integra Jarit instrument. Integra Jarit neither assumes nor authorizes any representative or other person to assume for it any other liability in connection with Integra Jarit's products.

All products are subject to slight changes in dimensions and finishes due to changes in manufacturing standards.

Ordering Procedure

A written purchase order on customer's form may be requested for purchase of Integra products.

Acceptance of Order

Orders are accepted upon approval by Integra.

Returned Goods Policy

Products must be returned in unopened packages with manufacturer's seals intact to be accepted for replacement or credit unless returned due to a complaint of product defect. Determination of a product defect will be made by Integra. Products will not be accepted for replacement if they have been in the possession of the customer for more than 90 days.

Repairs and Maintenance

Should your instruments require repair or maintenance, contact Integra for return authorization and address. Instruments returned to Integra for repair must have a statement testifying that each instrument has been thoroughly cleaned and sterilized. Failure to supply evidence of cleaning and disinfection will result in a cleaning charge and delayed processing of your instrument repair.

Product Information Disclosure

Integra and its subsidiaries ("Integra") and manufacturer exclude all warranties, except Integra's applicable standard warranty whether expressed or implied, including but not limited to, any implied warranties of merchantability or fitness for a particular purpose. Neither Integra nor manufacturer shall be liable for any incidental or consequential loss, damage, or expense, directly or indirectly arising from use of these products. Neither Integra nor manufacturer assume nor authorize any person to assume for them any other or additional liability or responsibility in connection with these products.

Disclaimer

All products are subject to slight changes in dimensions and finishes due to changes in manufacturing standards.

Care and Handling/IFU

For Care and Handling information and Instructions for use, please visit: <http://www.integralife.com/Jarit/Product-instructions.aspx?productID=589>

Limit uncertainty with exclusive tip patterns and flat suction instruments


Single Pivot Point

Single pivot video assisted thoracoscopic surgery instruments, like standard open surgical instruments, have one hinge. Single pivot instruments offer the best of user-friendly standard open instruments and the unique maneuverability of modern minimally invasive instrumentation.



Dual Pivot Point

Dual pivot video assisted thoracoscopic surgery instruments have two hinges with a shaft between. Unlike single pivot point instruments, dual pivot points offer unique access and functionality. The entire working length of the instrument, from the tip until the second pivot point, can be inserted through an incision or port while maintaining full maneuverability.



Tip Configurations



0.8mm Horizontal Serrations



1.2mm Horizontal Serrations



Vertical Serrations



Cross Serrations

Foerster Clamps — Dual Pivot

	Diameter	Length	0.8mm Horizontal Serrations	1.2mm Horizontal Serrations	Vertical Serrations	Cross Serrations
15° Bend	7mm Diameter	13¼" (335mm)	323-110	323-110A	323-110B	323-110C
	10mm Diameter	13¼" (339mm)	323-112	323-112A	323-112B	323-112C




7mm Diameter, 15° Bend




10mm Diameter, 15° Bend

Foerster Clamps — Dual Pivot



	Diameter	Length	0.8mm Horizontal Serrations	1.2mm Horizontal Serrations	Vertical Serrations	Cross Serrations
25° Bend	7mm Diameter	13" (332mm)	323-114	323-114A	323-114B	323-114C
	10mm Diameter	13¼" (336mm)	323-116	323-116A	323-116B	323-116C



Tip Configurations



11mm Triangular Ring



24mm Triangular Ring



Harken #2, 1x2 Teeth,
S-Curved (Open)



Harken #2, 1x2 Teeth,
S-Curved (Closed)

Duval Clamps — Dual Pivot

Catalog #	Description
323-120	13½" (344mm), 10mm, 11mm Triangular Ring, Curved Left
323-122	13¾" (348mm), 10mm, 24mm Triangular Ring, Curved Left



Duval Clamp— 11mm Triangular Ring



Duval Clamp— 24mm Triangular Ring

Cardio-Grip™ Harken Clamps — Dual Pivot

Catalog #	Description
323-150	Harken #1, 14" (353mm), 10mm, DeBakey Tapered Jaws, 1x2 Teeth, Slight Curve, No Ratchet
323-151	Harken #2, 13¾" (350mm), 10mm, DeBakey Tapered Jaws, 1x2 Teeth, S-Curved, No Ratchet



Cardio-Grip Harken #1 Clamp— Slight Curve, No Ratchet


Tip Configurations



1x2 Teeth, Curved Left



2x3 Teeth



1x2 Teeth, Curved Left



1x2 Teeth, Straight

Cardio-Grip™ Short DeBakey Clamps — Dual Pivot

Catalog #	Description
323-130	12¾" (325mm), 10mm, 1x2 Teeth, Curved Left
323-132	14" (355mm), 10mm, 1x2 Teeth, Curved Left
323-134	12¼" (310mm), 5mm, 1x2 Teeth, Curved Left
323-135	14" (355mm), 10mm, 2x3 Teeth, Curved Left



Cardio-Grip Short DeBakey Clamp— 10mm

Cardio-Grip Long DeBakey Clamps — Dual Pivot

Catalog #	Description
323-140	14½" (368mm), 10mm, Tapered Jaws, 1x2 Teeth, Curved Left
323-141	14¾" (373mm), 10mm, Tapered Jaws, 1x2 Teeth, Straight



Cardio-Grip Long DeBakey Clamp — 10mm

Tip Configurations



90° Angled



1.4mm DeBakey Jaws



2mm DeBakey Jaws

DeBakey Cooley Forceps — Dual Pivot

Catalog #	Description
323-202	13" (331mm), 10mm, 90° Angled, Ratchet
323-203	13" (331mm), 10mm, 90° Angled, No Ratchet



DeBakey Cooley Forceps — 10mm

Cardio-Grip™ DeBakey Forceps, Castroviejo Handle — Dual Pivot

Catalog #	Description
323-420	16½" (420mm), Castroviejo Handle, 5mm, 1.4mm DeBakey Jaws, No Ratchet
323-422	16¼" (411mm), Castroviejo Handle, 5mm, 2mm DeBakey Jaws, No Ratchet



Cardio-Grip DeBakey Forceps, Castroviejo Handle — 5mm

Tip Pattern



Short-Jaw Ring Forceps — Single Pivot

Catalog #	Description
323-015	12 1/2" (320mm), Straight
323-011	12 1/4" (308mm), Curved Jaw, Straight Handle
323-012	10 3/4" (275mm), "C" Shaped



Straight



Curved Jaw, Straight Handle




"C" Shaped

Tip Pattern



Short-Jaw Ring Forceps — Single Pivot

Catalog #	Description
323-013	12" (305mm), Curved On Side
323-014	12" (302mm), Curved Handle



Tip Pattern




Long-Jaw Ring Forceps — Single Pivot

Catalog #	Description
323-021	12 1/4" (310mm), Curved Jaw
323-022	11" (277mm), C-Shaped
323-023	12 1/4" (308mm), Curved on Side
323-024	12 1/4" (308mm), Curved Handle



Curved Jaw



C-Shaped



Curved on Side



Curved Handle

Tip Configurations



90° Ring




45° Ring



Straight

Long-Jaw Ring Forceps — Single Pivot

Catalog #	Description
323-025	12" (302mm), 90° Ring
323-026	12¼" (310mm), 45° Ring
323-027	12½" (320mm), Straight



90° Ring



45° Ring



Straight

Tip Pattern



Cardio-Grip™ Atraumatic Thoracoscopic Ring Forceps — Single Pivot

Catalog #	Description
323-053	12 3/4" (322mm), Straight
323-051	11 3/4" (298mm), "S" Shaped
323-052	12 1/4" (312mm), Curved



Straight



"S" Shaped



Curved

Tip Configurations



Slight Curve



Straight

Cardio-Grip™ Atraumatic Thoracoscopic Forceps — Single Pivot

Catalog #	Description
323-060	12" (304mm), Straight, No Ratchet
323-061	11¾" (300mm), Slight Curve, No Ratchet
323-062	12" (304mm), Straight, Ratchet
323-063	11¾" (300mm), Slight Curve, Ratchet



Straight, No Ratchet



Slight Curve, No Ratchet



Straight, Ratchet



Slight Curve, Ratchet

Tip Configurations



Cardio-Grip™ Atraumatic Allis Ring Forceps — Single Pivot

Catalog #	Description
323-070	12" (305mm), 2x3 Teeth, Straight



Cardio-Grip Atraumatic Allis Ring Forceps — Straight

Allis Forceps — Dual Pivot

Catalog #	Description
323-270	13½" (342mm), 7mm, 6x7 Teeth, Curved Left



Allis Forceps — 7mm, Curved Left

Tip Configurations



Curved



Straight



Curved



Straight

Carb-Bite™ Needle Holder — Dual Pivot

Catalog #	Description
323-640	13¼" (336mm), 7mm, Tungsten Carbide, Curved, Ratchet, For 4-o Suture
323-641	13¼" (336mm), 7mm, Tungsten Carbide, Straight, Ratchet, For 4-o Suture



7mm, for 4-o Suture

Carb-Bite Needle Holder — Dual Pivot

Catalog #	Description
323-620	13¼" (336mm), 10mm, Tungsten Carbide, Curved, Ratchet, For 2-o Suture
323-621	13¼" (336mm), 10mm, Tungsten Carbide, Straight, Ratchet, For 2-o Suture



10mm, for 2-o Suture

Tip Configurations



Straight




Curved Left



Straight

Jarit® Stardust® Needle Holder — Dual Pivot

Catalog #	Description
323-651	11" (277mm), 5mm, Jarit Stardust, Straight, Ratchet, For 5-o Suture



Jarit Stardust Needle Holder — 5mm, For 5-o Suture

Jarit Stardust Castroviejo Needle Holder — Dual Pivot

Catalog #	Description
323-650	15½" (395mm), 5mm, Jarit Stardust, Curved Left, Ratchet, For 5-o Suture
323-653	15¼" (387mm), 5mm, Jarit Stardust, Straight, Ratchet, For 5-o Suture



Jarit Stardust Castroviejo Needle Holder — 5mm, Curved Left, For 5-o Suture



Jarit Stardust Castroviejo Needle Holder — 5mm, Straight, For 5-o Suture

Tip Configurations



Blunt Micro Blades, Straight



Blunt Blades, Curved



Blunt Blades, Curved

Jarit® Supercut™ Scissors — Dual Pivot

Catalog #	Description
323-701	9" (227mm), 5mm, Blunt Micro Blades, Straight

Jarit Supercut Scissors — 5mm

Jarit Supercut Dissecting Scissors — Dual Pivot




Catalog #	Description
323-720	11¼" (287mm), 5mm, Blunt Blades, Curved

Jarit Supercut Dissecting Scissors — 5mm

Jarit Supercut Metzenbaum Scissors — Dual Pivot

Catalog #	Description
323-730	11¾" (296mm), 10mm, Blunt Blades, Curved

Jarit Supercut Metzenbaum Scissors — 10mm




Tip Pattern



Curved Left

Dennis Dissector — Dual Pivot

Catalog #	Description
323-280	13" (328mm), 7mm, Serrated, Curved Left



Dennis Dissector — 7mm

Tip Configurations



Oblong Ring



Open Cup

Node Graspers, Oblong Ring — Dual Pivot

Catalog #	Description
323-100	13¼" (338mm), 5mm, Oblong Ring Jaw, Curved Left
323-101	13" (330mm), 5mm, Oblong Ring Jaw, Straight



Node Grasper, Curved Left, Oblong Ring — 5mm



Node Grasper, Straight, Oblong Ring — 5mm

Node Graspers, Open Cup — Dual Pivot

Catalog #	Description
323-102	13½" (341mm), 10mm, Open Cup Jaw, Curved Left
323-103	13½" (344mm), 10mm, Open Cup Jaw, Straight



Node Grasper, Curved Left, Open Cup — 10mm



Node Grasper, Straight, Open Cup — 10mm

Tip Configurations



5mm, Straight



5mm, Curved, Flat



10mm, Straight



10mm, Curved, Flat

Suction Instruments

Catalog #	Description
323-900	18 ³ / ₄ " (475mm), 5mm, Blunt Tip, Straight
323-901	18 ³ / ₄ " (475mm), 5mm, Variable Suction, Blunt Tip, Straight
323-902	16 ¹ / ₄ " (413mm), 5mm, Flat Suction, Blunt Tip, Curved
323-903	16 ¹ / ₄ " (413mm), 5mm, Variable Flat Suction, Blunt Tip, Curved
323-910	18 ³ / ₄ " (475mm), 10mm, Blunt Tip, Straight
323-911	18 ³ / ₄ " (475mm), 10mm, Variable Suction, Blunt Tip, Straight
323-912	16 ¹ / ₂ " (421mm), 10mm, Flat Suction, Blunt Tip, Curved
323-913	16 ¹ / ₂ " (421mm), 10mm, Variable Flat Suction, Blunt Tip, Curved




Blunt Tip Variable Suction — 5mm, Straight



Blunt Tip Variable Flat Suction — 5mm, Curved



Blunt Tip Variable Suction — 10mm, Straight



Blunt Tip Variable Flat Suction — 10mm, Curved

Tip Configurations




Curved




Curved

Multi-Hole Tip Suction Instruments

Catalog #	Description
323-920	12" (307mm), 3.2mm, Multi-Hole Tip, Curved
323-921	12" (307mm), 3.2mm, Variable Suction, Multi-Hole Tip, Curved
323-922	14" (357mm), 3.2mm, Multi-Hole Tip, Curved
323-923	14" (357mm), 3.2mm, Curved Variable Suction, Multi-Hole Tip




Multi-Hole Tip, Curved, 14"



Multi-Hole Tip Variable Suction, Curved, 12"

Dennis Suction Instruments

Catalog #	Description
323-914	14½" (369mm), 5mm, Multi-Hole Tip, Curved
323-915	14½" (369mm), 5mm, Variable Suction, Multi-Hole Tip, Curved
323-916	18½" (469mm), 5mm, Multi-Hole Tip, Curved
323-917	18½" (469mm), 5mm, Variable Suction, Multi-Hole Tip, Curved



Variable Dennis Suction—Multi-Hole Tip, Curved, 18½"



Variable Dennis Suction—Multi-Hole Tip, Curved, 14½"

Tip Configurations




Curved



Angled

Yankauer Suction Instruments


Catalog #	Description
323-930	13 ¹ / ₄ " (335mm), Curved, With Thumb Valve
323-TIP-930	Tip for Yankauer Suction Tube, 1/pkg



Yankauer Suction — Curved, with Thumb Valve

Andrews Suction Instruments

Catalog #	Description
323-931	11 ³ / ₄ " (300mm), Angled
323-TIP-931	Tip for Andrews Suction Tube, 1/pkg



Andrews Suction — Angled

Brush Set for Suction Instruments

Catalog #	Description
323-990	19" (480mm), 5mm, For Straight Suction, Pack of 10
323-991	19 ³ / ₄ " (500mm), 10mm, For Straight Suction, Pack of 10



Tip Configurations



Knot Tier/Pusher



Supercut Knot Tier/Pusher

Knot Tier/Pushers — Dual Pivot


Catalog #	Description
323-800	11" (280mm), 5mm, For 4-o Suture
323-801	12" (305mm), 5mm, For 4-o Suture
323-810	12½" (318mm), 5mm, For 2-o Suture
323-811	15¼" (390mm), 5mm, For 2-o Suture



Knot Tier/Pushers — 5mm

Supercut Knot Tier/Pusher — Dual Pivot

Catalog #	Description
323-820	12¼" (310mm), 5mm, For 2-o Suture



Supercut Knot Tier/Pusher — 5mm



Limit Uncertainty

with the Brands you Trust.

323-011.....	9	323-112B.....	4	323-650.....	17	323-TIP-931.....	23
323-012.....	9	323-112C.....	4	323-651.....	17	323-990.....	23
323-013.....	10	323-114.....	5	323-653.....	17	323-991.....	23
323-014.....	10	323-114A.....	5	323-701.....	18		
323-015.....	91	323-114B.....	5	323-720.....	18		
323-021.....	11	323-114C.....	5	323-730.....	18		
323-022.....	11	323-116.....	5	323-800.....	24		
323-023.....	11	323-116A.....	5	323-801.....	24		
323-024.....	11	323-116B.....	5	323-810.....	24		
323-025.....	12	323-116C.....	5	323-811.....	24		
323-026.....	12	323-120.....	6	323-820.....	24		
323-027.....	12	323-122.....	6	323-900.....	21		
323-051.....	13	323-130.....	7	323-901.....	21		
323-052.....	13	323-132.....	7	323-902.....	21		
323-053.....	13	323-134.....	7	323-903.....	21		
323-060.....	14	323-135.....	7	323-910.....	21		
323-061.....	14	323-140.....	7	323-911.....	21		
323-062.....	14	323-141.....	7	323-912.....	21		
323-063.....	14	323-150.....	6	323-913.....	21		
323-070.....	15	323-151.....	6	323-914.....	22		
323-100.....	20	323-202.....	8	323-915.....	22		
323-101.....	20	323-203.....	8	323-916.....	22		
323-102.....	20	323-270.....	15	323-917.....	22		
323-103.....	20	323-280.....	19	323-920.....	22		
323-110.....	4	323-420.....	8	323-921.....	22		
323-110A.....	4	323-422.....	8	323-922.....	22		
323-110B.....	4	323-620.....	16	323-923.....	22		
323-110C.....	4	323-621.....	16	323-930.....	23		
323-112.....	4	323-640.....	16	323-TIP-930.....	23		
323-112A.....	4	323-641.....	16	323-931.....	23		

A

Allis Forceps	15
Andrews Suction	23
Atraumatic Allis Ring Forceps	15
Atraumatic Thoracoscopic Forceps	14
Atraumatic Thoracoscopic Ring Forceps	13

B

Brush Set	23
-----------------	----

C

Carb-Bite Needle Holder	16
Castroviejo Needle Holder	17
Clamps	4-7

D

DeBakey Clamps	7
DeBakey Cooley Forceps	8
DeBakey Forceps	8
Dennis Dissector	19
Dennis Suction	22
Dissecting Scissors	18
Dissector	19
Duval Clamps	6

F

Flat Suction	21
Foerster Clamps	4-5
Forceps	8-10, 13, 15

H

Harken Clamps	6
---------------------	---

K

Knot Tier/Pusher	24
------------------------	----

L

Long DeBakey Clamps	7
Long-Jaw Ring Forceps	11-12

M

Metzenbaum Scissors	18
Multi-Hole Tip Suction	22

N

Needle Holders	16-17
Node Grasper	20

R

Ring Forceps	9-15
--------------------	------

S

Scissors	18
Short DeBakey Clamps	7
Short-Jaw Ring Forceps	8-10
Stardust Castroviejo Needle Holder	17
Stardust Needle Holder	17
Suction	21-23
Supercut Dissecting Scissors	18
Supercut Knot Tier/Pusher	24
Supercut Metzenbaum Scissors	18
Supercut Scissors	18

T

Thoracoscopic Forceps	14
-----------------------------	----

Y

Yankauer Suction	23
------------------------	----

For more information or to place an order, please contact:

877-468-5572 USA ■ 800-275-7192 fax

integralife.com



Integra, the Integra logo, Jarit and Jarit Stardust are registered trademarks of Integra LifeSciences Corporation or its subsidiaries in the United States and/or other countries. Jarit Supercut, Cardio-Grip, Carb-Bite and Carb-Edge are trademarks of Integra LifeSciences Corporation or its subsidiaries.

All other trademarks are the property of their respective owners. ©2013 Integra LifeSciences Corporation. All rights reserved.

Printed in the USA. 2.5M 0075179-2