


Integra®

Ruggles®-Redmond™

Spinal Curettes

Limit uncertainty with spinal
curettes in the sizes, coatings,
and configurations you need



INTEGRA®
LIMIT UNCERTAINTY

Ultra-Coat® Spinal Fusion Curettes


Ultra-Coat Spinal Fusion Curettes have a hard material coating that increases the surface hardness versus regular stainless steel curettes.



Cup Size	Straight				Angled			Reverse Angled		
	6¾" (171mm)	8" (203mm)	9" (229mm)	12" (305mm)	8" (203mm)	9" (229mm)	12" (305mm)	8" (203mm)	9" (229mm)	12" (305mm)
5-0	RB7701	RB7712	RB7723	RB7734	RB7745	RB7756	RB7767	RB7801	RB7812	RB7823
4-0	RB7702	RB7713	RB7724	RB7735	RB7746	RB7757	RB7768	RB7802	RB7813	RB7824
3-0	RB7703	RB7714	RB7725	RB7736	RB7747	RB7758	RB7769	RB7803	RB7814	RB7825
2-0	RB7704	RB7715	RB7726	RB7737	RB7748	RB7759	RB7770	RB7804	RB7815	RB7826
0	RB7705	RB7716	RB7727	RB7738	RB7749	RB7760	RB7771	RB7805	RB7816	RB7827
#1	RB7706	RB7717	RB7728	RB7739	RB7750	RB7761	RB7772	RB7806	RB7817	RB7828
#2	RB7707	RB7718	RB7729	RB7740	RB7751	RB7762	RB7773	RB7807	RB7818	RB7829
#3	RB7708	RB7719	RB7730	RB7741	RB7752	RB7763	RB7774	RB7808	RB7819	RB7830
#4	RB7709	RB7720	RB7731	RB7742	RB7753	RB7764	RB7775	RB7809	RB7820	RB7831
#5	RB7710	RB7721	RB7732	RB7743	RB7754	RB7765	RB7776	RB7810	RB7821	RB7832
#6	RB7711	RB7722	RB7733	RB7744	RB7755	RB7766	RB7777	RB7811	RB7822	RB7833

Apfelbaum Extended Reach Curettes

These curettes offer an extended reach to more easily remove posterior osteophytes. Each reach length (5mm, 6mm, and 7mm) features a color-coded handle for easy identification.



Cup Size	Reach	12" (305mm)	14½" (368mm)
2mm x 3mm	5mm	RN1171	RN1181
2mm x 3mm	6mm	RN1172	RN1182
2mm x 3mm	7mm	RN1173	RN1183

Cervical and Lumbar Bayonet Curettes

The bayonet handle may help improve visualization.



Cervical Bayonet Curettes - 8½"

Cup Size	Forward Straight	Forward Angle	Backward Straight	Backward Angle
4-O	RR8113	RR8118	RR8123	RR8128
3-O	RR8112	RR8117	RR8122	RR8127
2-O	RR8111	RR8116	RR8121	RR8126
1-O	RR8110	RR8115	RR8120	RR8125
2-O BC	RR8114	RR8119	RR8124	RR8129

Catalog #	Description	Qty
RR8100	Cervical Bayonet Curette Set 8½"	
RR8130	Cervical Bayo Nerve Hook 8¾"	1
RR8131	Cervical Ultra Bayo Nerve Hook 8¾"	1
RR8132	Cervical Bayo Waxer 8¾"	1
RR8133	Cervical Bayo Knife Handle	1
RS0900	Bayonet Curette Tray	2

XL Cervical Bayonet Curettes - 9½"

Cup Size	Forward Straight	Forward Angle	Backward Straight	Backward Angle
4-O	RR8213	RR8218	RR8223	RR8228
3-O	RR8212	RR8217	RR8222	RR8227
2-O	RR8211	RR8216	RR8221	RR8226
1-O	RR8210	RR8215	RR8220	RR8225
2-O BC	RR8214	RR8219	RR8224	RR8229


Catalog #	Description	Qty
RR8200	XL Cervical Bayonet Curette Set 9½"	
RR8230	XL Cervical Bayo Nerve Hook 9½"	1
RR8231	XL Cervical Ultra Bayo Nerve Hook 9½"	1
RR8232	XL Cervical Waxer 9½"	1
RR8133	Cervical Bayo Knife Handle	1
RS0900	Bayonet Curette Tray	2

Lumbar Bayonet Curettes - 9½"

Cup Size	Forward Straight	Forward Angle	FD	Backward Straight	Backward Angle	BD
1-O	RR8010	RR8013	RR8016	RR8019	RR8022	RR8025
#1	RR8011	RR8014	RR8017	RR8020	RR8023	RR8026
#2	RR8012	RR8015	RR8018	RR8021	RR8024	RR8027

Catalog #	Description	Qty
RR8000	Lumbar Bayonet Curette Set 9½"	
RR8028	Lumbar Bayo Penfield Dissector #1	1
RR8029	Lumbar Bayo Penfield Dissector #2	1
RR8030	Lumbar Bayo Nerve Hook 9½"	1
RR8031	Lumbar Bayo Penfield Dissector #4	1
RR8032	Lumbar Bayo Waxer 9½"	1
RR8033	Lumbar Bayo Spine Probe 9½"	1
RS0900	Bayonet Curette Tray	2


Curette Cup Sizes *Size chart is for visual approximation, not actual size.




Flatback Spinal Curettes

The flatback cup's slimmer profile design helps provide improved maneuverability versus standard cup design when removing disc material.

Flatback Spinal Curette



Standard Spinal Curette



Standard

Cup Size	Straight	45° (9") Angled	90° (8¾") Angled	Reverse Angled
6-O	R2799	R2819	R2839	RB7170
5-O	R2800	R2820	R2840	R2850
4-O	R2801	R2821	R2841	R2851
3-O	R2802	R2822	R2842	R2852
2-O	R2803	R2823	R2843	R2853
O	R2804	R2824	R2844	R2854
#1	R2805	R2825	R2845	R2855
#2	R2806	R2826	R2846	R2856
#3	R2807	R2827	R2847	R2857
#4	R2808	R2828	R2848	R2858
#5	R2809	R2829	R2849	R2859
#6	R2810	R2830	R2850	RB7181

Ultra-Coat®

Straight	45° (9") Angled	90° (8¾") Angled	Reverse Angled
RB7900	RB7920	RB7940	RB7960
RB7901	RB7921	RB7941	RB7961
RB7902	RB7922	RB7942	RB7962
RB7903	RB7923	RB7943	RB7963
RB7904	RB7924	RB7944	RB7964
RB7905	RB7925	RB7945	RB7965
RB7906	RB7926	RB7946	RB7966
RB7907	RB7927	RB7947	RB7967
RB7908	RB7928	RB7948	RB7968
RB7909	RB7929	RB7949	RB7969
RB7910	RB7930	RB7950	RB7970
RB7911	RB7931	RB7951	RB7971

Other Ruggles-Redmond Spinal Curettes

- Cone Ring Curettes
- Cobb Style Curettes
- Scoville Disc Curettes
- Epstein Curettes
- Semmes Spinal Curettes
- Sisco Bayonet Curettes
- MicroAngled Curettes
- MicroBayonet Curettes
- Round Handle Bayonet Micro Curettes
- ThoracoLumbar Curettes, Resin Handles

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Please consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions, and instructions for use.
- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.

For more information or to place an order, please contact:

United States, Canada, Asia, Pacific, Latin America

USA 800-654-2873 ▪ 888-980-7742 fax

International +1 609-936-5400 ▪ +1 609-750-4259 fax

integralife.com/contact



Integra York PA, Inc.
589 Davies Drive
York, PA 17402 ▪ USA