

# RF15A



Portable High-Definition Digital Bronchoscope




## Versatility

It is suitable for trachea and difficult airway intubation, airway and double-lumen inspection. Alveolar lavage, biopsy sampling, pneumothorax closure, removal of tracheal foreign bodies and other bronchial treatments.

## Applicability

Application departments: ICU, NICU, respiratory department, anesthesiology department, emergency department, etc., basically meet the clinical needs of each department



### Optical system

■ Field of view	90°
■ Depth of field	3–50mm
■ Direction of view	0°

### Insertion portion

■ Insertion tube OD	4.0mm/4.9mm/5.5mm
■ Angulation range	U 180° / D 130°
■ Inside diameter	1.2mm/2.2mm/2.6mm
■ Working length	600mm
■ Max. resolution	1280×720P



**ZHONGSHAN WESEE MEDITECH CO., LTD.**

Address: 6F, Block C, 2Jianshuo Str. HoChung Village, Nanlang, Zhongshan 528451 P. R. China

TEL: (+86)0760-86960005

Email: info@weseemed.com

www.weseemed.com

